

Connectivity options for Audio Service hearing aids

Form factor			Audio streaming				Hands-free calling			Control via Audio Service app
			Bluetooth® Classic	Bluetooth® LE Audio	ASHA	MFi	Bluetooth® Classic	Bluetooth® LE Audio	MFi	
G8	RIC	R Li T BC	●	-	-	-	●	-	-	●
	RIC	R Li, R Li T	-	●	●	●	-	●	●	●
	SR	SR Li	-	●	●	●	-	●	●	●
	BTE	B M, B P Li, B SP Li	-	●	●	●	-	●	●	●
	IF	IF CIC Li	-	-	-	-	-	-	-	●
	Custom	CIC Li	-	-	-	-	-	-	-	●
	Custom	DIC, CIC, ITC, ITE	-	-	-	-	-	-	-	●
G7	RIC	R S, R Li, R Li T	-	-	●	●	-	-	●	●
	SR	SR Li	-	-	●	●	-	-	●	●
	BTE	B M, B P, B SP	-	-	●	●	-	-	-	●
	Custom	ITC Li, ITE Li	-	-	●	●	-	-	●	●
G4	RIC	R Li	-	-	-	-	-	-	-	●
	BTE	B P Li, B SP Li	-	-	-	-	-	-	-	●

Product availability and compatibility may vary depending on country and Tech Level.
 The Bluetooth® word mark and logos are the property of Bluetooth SIG Inc.
 Any use of such marks by WSAUD A/S is governed by a license agreement.
 Other trademarks and trade names remain the property of their respective owners.